

Comprefit®

Strongest Containment



Compreflex®

Easy & Measurable



Coolflex™

Most Breathable

Inelastic Wraps Guide

LOWER EXTREMITY



Simple Solutions for

Venous Disease, Lymphedema & Wound Care

Compreflex Standard

CALF & FOOT







When to Use

- · Adjustable ranges of resting compression is desired (20-50mmHg)
- · Non-attached foot piece is desired
- · No compression in the liner is desired
- · Patient requires a calf and foot garment

Ideal for

- Chronic Venous Disease¹
- · Chronic Venous Insufficiency (CVI) including, but
- · not limited to: ulcer due to venous stasis
- · Lymphedema Active and Maintenance
- · Non-specific or post-surgical edema



Compreboot® Standard

Variable compression



Basic Liners Non-Compressive



Black

Beige

Accutab® Devices

Compreflex Standard

CALF

20-50mmHg





When to Use

- · Patient requires compression in the calf only or a different foot option
- · Varying ranges of adjustable resting compression is desired (20-50mmHg)
- · No compression in the liner is desired
- · No Compreboot needed
- · Ability to wear standard footwear over liner is desired

Ideal for

- · Lymphedema
- · Edema
- · Chronic venous disease (CVD)
- · Chronic venous insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- · Lipedema













Wrap Details

20-50mmHg Compression

Each Compreflex Calf wrap includes the Accutab, allowing users to accurately adjust the garment to the desired level of compression

Intelligent Design

Garment features a contoured calf, providing a superior fit

Easy Application

With a front stretch panel to secure the garment in place and straps that roll back, the garment is easy to don and adjust, promoting adherence

Premium Materials

Soft, conforming, breathable material provides a comfortable fit

Superior Quality Each garment is double stitched to ensure durability and is backed by a six-month warranty**



Comprefit Standard

CALF & FOOT

20-50mmHg





Garment Features

- · 20-50mmHg firm, non-stretch compression
- Rigid, non-stretch garment profile providing a high level of containment
- · Two slot options provide a comfortable fit
- Garment can be adjusted to include an additional 10cm of circumferential range by adding the Comprefit Standard Strap Extenders

Ideal for

- · Chronic Venous Disease
- Chronic Venous Insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- · Lymphedema
- · Non-specific or post surgical edema

When to Use

- · Patient requires a garment with no stretch
- Patient requires compression in foot; product includes Compreboot
- · A high level of containment is needed
- · No compression in the liner is desired





Pair of Basic Liners

Non-Compressive





Black







Coolflex Standard

CALF

20-50mmHa





Garment Features

- · 20-50mmHg of lightweight, gradient compression
- Garment is made of Spacer Fabric and includes superior properties preventing the trapping of heat and promoting the transfer of moisture
- EZ Close Pressure™ system visually notes the compression levels of 20–30, 30–40, and 30–50mmHg on the garment straps for the easy application of compression
- · Hook closures for maximum strength and durability
- One-handed side closure system promoting independence and patient adherence

Ideal for

- · Chronic Venous Disease
- Chronic Venous Insufficiency (CVI) including, but not limited to: ulcer due to venous stasis
- Lymphedema
- · Non-specific or post surgical edema

When to Use

- · Patient requires a lightweight material due to skin sensitivities or fragility
- Varying ranges of adjustable resting compression is desired (20–50mmHg)
- · Patient has limited hand strength or dexterity
- Ability to wear standard footwear over liner is desired
- Patient has difficulty leaning forward or bending, lateral or side closure to aid the patient in applying the garment is needed









Wrap Details

20-50mmHg Compression

Coolflex is a lightweight compression device optimal for individuals with fragile or broken skin. Its patented EZ-Close pressure system closes laterally which can promotes independence and improves patient outcomes.

Most Breathable Garment

The Spacer Fabric is designed to be highly breathable for cool-wearing comfort

One-Handed Donning

EZ-Close Pressure System of eight tabs makes one-handed, lateral donning simple, even for patients with dexterity issues


Innovative Design

Durable, metal hook closures for maximum strength and durability



The integrated EZ-Close Pressure System provides a visual indicators for 20–30, 30–40, and 30–50mmHg on each tab.





Please consult your doctor, or other qualified health care professional before using any product discussed within the literature. See product insert for complete information.

* Garment includes a 6-month warranty against defects in materials and workmanship