

Graduated compression
socks and hosiery

SIGVARIS
GROUP



SIGVARIS GROUP WELL BEING

Care for Yourself



Choice is empowering. That's why SIGVARIS GROUP offers styles and patterns that let you be you. We use the world's most premium yarns so you can express your personality while enjoying the clinically proven benefits of compression every day, no matter what the day brings.

**With over 150 years of success,
we are looking towards the future.**

We are renowned for making exceptional products that combine advanced technology like double covered yarns and quality craftsmanship with diverse styles that fit our modern customers' lives.

Our reason: **You.** We want you to live a better life, do what you love and feel your best. So, go further with every step in SIGVARIS GROUP graduated compression.



We are devoted
to giving you the **freedom,
confidence, and security**
to live your life to the fullest.

SIGVARIS GROUP WELL BEING

Care for Yourself

Sigvaris Well Being is the easy way to get the energizing feeling of compression. Sigvaris Well Being socks and hosiery make shopping for compression simple with a sizing system based on your shoe size.

Ideal for:

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease, prevention of swelling and varicose veins
- Simplified sizing

Table of contents

Series / Product	Page
146 / 186 Casual Cotton	3
151 / 191 Sea Island Cotton	5
152 / 192 All-Season Merino Wool	7
143 / 183 Microfiber Shades	9
120 Sheer Fashion	11
170 Soft Silhouette Leggings	13
401 Athletic Recovery Socks	15
412 Performance Sleeves	17
421 Merino Outdoor Socks	19
160 Eversoft Diabetic Socks	21
Accessories	23
Dealer Portal	25
How to order	26



SIGVARIS
GROUP

CASUAL COTTON 146/186

for women and men

Compression Level

15–20mmHg


Features and Benefits

- » Made with premium Supima cotton for optimum softness and durability
- » Flat, low-friction toe seam to prevent irritation
- » Soft, stretchable, knit-in band prevents pinching

Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease

Color Swatches



Fiber Content

- » 66% Supima Cotton, 26% Nylon, 8% Spandex
- » Non-latex

SUPIMA®
WORLD'S FINEST COTTONS



146CB30
15–20mmHg
Closed Toe Calf Size B in Khaki



146CC11
15–20mmHg
Closed Toe Calf Size C in Brown



186CB00
15–20mmHg
Closed Toe Calf Size B in White



186CB99
15–20mmHg
Closed Toe Calf Size B in Black

146 Women's Calf Sizing Chart

Women's Shoe	5–7	7.5–9.5	10–12
Ankle Circumference	7"–9" (18–23cm)	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)
Size	A	B	C

Product Code Example: 146CB10

Series	Style	Size	Color
146	C	B	10

186 Men's Calf Sizing Chart

Men's Shoe	6.5–8.5	9–11	11.5–14
Ankle Circumference	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)	8.5"–11.5" (21–29cm)
Size	A	B	C

Product Code Example: 186CA30

Series	Style	Size	Color
186	C	A	30

Styles and Colors

Closed Toe = ●	White 00	Navy 10	Brown 11	Khaki 30	Black 99
Calf (C)	●	●	●	●	●



SIGVARIS
GROUP

SEA ISLAND COTTON 151/191

for women and men

Compression Level

15–20mmHg

Features and Benefits

- » Knitted from the world's rarest cotton
- » Sea Island Cotton yarns are lightweight and have great strength, a silky luster, and an incredibly soft feel
- » Breathable and comfortable for sensitive skin
- » Ideal for travel and everyday wear


Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease

Fiber Content

- » 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- » Non-latex

Color Swatches



SUPIMA®
WORLD'S FINEST COTTONS



151CA10
15–20mmHg
Closed Toe Calf Size A in Navy



191CB99
15–20mmHg
Closed Toe Calf Size B in Black

151 Women's Calf Sizing Chart

Women’s Shoe	5–7	7.5–9.5	10–12
Ankle Circumference	7”–9” (18–23cm)	7.5”–10” (19–25cm)	8”–10.5” (20–27cm)
Size	A	B	C

Product Code Example: 151CA10

Series	Style	Size	Color
151	C	A	10

191 Men's Calf Sizing Chart

Men’s Shoe	6.5–8.5	9–11	11.5–14
Ankle Circumference	7.5”–10” (19–25cm)	8”–10.5” (20–27cm)	8.5”–11.5” (21–29cm)
Size	A	B	C

Product Code Example: 191CB99

Series	Style	Size	Color
191	C	B	99

Styles and Colors

Closed Toe = ●	Navy 10	Black 99
Calf (C)	●	●

ALL-SEASON 152/192 MERINO WOOL

for women and men

Compression Level

15–20mmHg

Features and Benefits

- » Soft and comfortable
- » Natural properties of wool make these socks odor resistant
- » Thermoregulating for year-round wear and keeps feet dry


Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease

Fiber Content

- » 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- » Non-latex

Color Swatches





152CC99
15–20mmHg
Closed Toe Calf Size C in Black



192CB10
15–20mmHg
Closed Toe Calf Size B in Navy

152 Women's Calf Sizing Chart

Women’s Shoe	5–7	7.5–9.5	10–12
Ankle Circumference	7”–9” (18–23cm)	7.5”–10” (19–25cm)	8”–10.5” (20–27cm)
Size	A	B	C

Product Code Example: 152CA99

Series	Style	Size	Color
152	C	A	99

192 Men's Calf Sizing Chart

Men’s Shoe	6.5–8.5	9–11	11.5–14
Ankle Circumference	7.5”–10” (19–25cm)	8”–10.5” (20–27cm)	8.5”–11.5” (21–29cm)
Size	A	B	C

Product Code Example: 192CC10

Series	Style	Size	Color
192	C	C	10

Styles and Colors

Closed Toe = ●	Navy 10	Brown 11	Black 99
Calf (C)	●	●	●

MICROFIBER SHADES 143/183

for women and men

Compression Level

15–20mmHg

Features and Benefits

- » Advanced microfiber technology for ultimate comfort
- » Glides on and off with ease
- » Excellent breathability and moisture management
- » Soft and durable sock for everyday wear


Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease

Fiber Content

- » 87% Nylon, 13% Spandex
- » Non-latex

Color Swatches





143CC02
15–20mmHg
Closed Toe Calf Size C in
Mariner Stripe



183CB47
15–20mmHg
Closed Toe Calf Size B in
Graphite Argyle



143CB48
15–20mmHg
Closed Toe Calf Size B in Royal
Blue Argyle



183CB17
15–20mmHg
Closed Toe Calf Size B in
Graphite Stripe

143 Women's Calf Sizing Chart

Women's Shoe	5–7	7.5–9.5	10–12
Ankle Circumference	7"–9" (18–23cm)	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)
Size	A	B	C

Product Code Example: 143CB02

Series	Style	Size	Color
143	C	B	02

183 Men's Calf Sizing Chart

Men's Shoe	6.5–8.5	9–11	11.5–14
Ankle Circumference	7.5"–10" (19–25cm)	8"–10.5" (20–27cm)	8.5"–11.5" (21–29cm)
Size	A	B	C

Product Code Example: 183CC17

Series	Style	Size	Color
183	C	C	17

Styles and Colors

	Mariner Stripe 02	Pink Stripe 15	Graphite Stripe 17	Onyx Stripe 19	Graphite Heather 27	Purple Argyle 46	Graphite Argyle 47	Royal Blue Argyle 48
Calf (C)				●	●		●	●
Calf (C) <i>for women only</i>	●	●				●		
Calf (C) <i>for men only</i>			●					

SHEER FASHION 120

for women

Compression Level

15–20mmHg

Features and Benefits

- » Fabric feels soft against the skin
- » Incredible sheer design in a variety of fashionable colors
- » Sheer toe and heel area
- » Comfortable stretch is excellent for everyday wear
- » Cotton-lined gusset on pantyhose styles


Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease

Fiber Content

- » 84% Nylon, 16% Spandex (Calf style)
- » Content can vary by style
- » Non-latex

Color Swatches





120CC12
15–20mmHg
Closed Toe Calf Size C in Charcoal



120PC73
15–20mmHg
Closed Toe Pantyhose Size C in Café

Calf / Thigh Sizing Chart

Women’s Shoe	5–7	7.5–9.5	10–12
Ankle Circumference	7.5”–9” (19–23cm)	8”–10” (20.5–25cm)	8.5”–10.5” (22–27cm)
Size	A	B	C

Weight Guidelines For Thigh Style: Size A: <135Lbs., Size B: <180Lbs., Size C: <230Lbs.

Pantyhose Sizing Chart

	Lbs/Kg	100 / 45	105 / 48	110 / 50	115 / 52	120 / 54	125 / 57	130 / 59	135 / 61	140 / 64	145 / 66	150 / 68	155 / 70	160 / 73	165 / 75	170 / 77	175 / 80	180 / 82	185 / 84	190 / 86	195 / 89	200 / 91	205 / 93	210 / 95	215 / 98	220 / 100	225 / 102	230 / 105	
Ft/Cm																													
4'9" / 145cm	Size A Ankle Circumference 7.5"—9" (19–23cm)																												
4'10" / 147cm																													
4'11" / 150cm																													
5'0" / 152cm																													
5'1" / 155cm																													
5'2" / 157cm																													
5'3" / 160cm	Size B Ankle Circumference 7.5"—9" (19–23cm)																												
5'4" / 162cm																													
5'5" / 165cm																													
5'6" / 167cm																													
5'7" / 170cm																													
5'8" / 173cm																													
5'9" / 175cm	Size C Ankle Circumference 8"—10" (20.5–25cm)																												
5'10" / 178cm																													
5'11" / 180cm																													
6' / 183cm																													
	Size D Ankle Circumference 8"—10" (20.5–25cm)																												
	Size E Ankle Circumference 8.5"—10.5" (22–27cm)																												
	Size F Ankle Circumference 8.5"—10.5" (22–27cm)																												

Styles and Colors

Closed Toe = ●
Open Toe = ○

	Charcoal 12	Taupe 29	Honey 33	Golden 36	Café 73	Black 99
Calf (C)	●	●○	●○	●○	●	●
Thigh (N)	●	●	●	●	●	●
Pantyhose (P)	●	●	●	●	●	●

Product Code Example: 120CC036

Series	Style	Size	Open toe	Color
120	C	C	O	36

Product Code Example: 120NC99

Series	Style	Size	Color
120	N	C	99

SOFT SILHOUETTE LEGGINGS 170

for women

Compression Level

15–20mmHg

Features and Benefits

- » Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts
- » Foldable, comfort-stretch waistband provides a seamless look for all-day wearing comfort
- » Advanced soft stretch microfiber fabric shapes legs and smooths skin


Best Suited for

- » Tired, heavy, restless, aching legs
- » Travelers, occupational leg health risks, expecting mothers, family history of chronic venous disease

Fiber Content

- » 83% Nylon, 17% Spandex
- » Non-latex

Color Swatches



*Made on demand colors may take up to 7 days to receive



170LB99
15–20mmHg
Legging Size B in Black



Non-binding ankle cuff
stays comfortably in place
for all day wear

Leggings Sizing Chart

	Ft/Cm	Lbs/Kg	100 / 45	105 / 48	110 / 50	115 / 52	120 / 54	125 / 57	130 / 59	135 / 61	140 / 64	145 / 66	150 / 68	155 / 70	160 / 73	165 / 75	170 / 77	175 / 80	180 / 82	185 / 84	190 / 86	195 / 89	200 / 91	205 / 93	210 / 95	215 / 98	220 / 100	225 / 102	230 / 105
	4'10" / 147cm																												
	4'11" / 150cm																												
	5'0" / 152cm																												
	5'1" / 155cm																												
	5'2" / 157cm																												
	5'3" / 160cm																												
	5'4" / 162cm																												
	5'5" / 165cm																												
	5'6" / 167cm																												
	5'7" / 170cm																												
	5'8" / 173cm																												
	5'9" / 175cm																												
	5'10" / 178cm																												
	5'11" / 180cm																												
	6' / 183cm																												

Guidelines for Ankle Circumference
Size A: 7"–9" (18–23cm), Size B & C: 8"–10" (20–26cm), Size D: 8.5"–11.5" (22–29cm)

Styles and Colors

Open Foot =	Midnight Blue 09	Mulberry 60	Espresso 89	Black 99
Leggings (L)				

Product Code Example: 170LB99

Series	Style	Size	Color
170	L	B	99

ATHLETIC 401 RECOVERY SOCKS

for women and men

Compression Level

15–20mmHg

Features and Benefits

- » Helps reduce lactic acid buildup¹
- » Comfortable and easy to put on and take off
- » *drirelease*® yarns provide moisture wicking and odor control
- » No compression in heel and instep for ease of donning


Best Suited for

- » Tired, heavy, restless, aching legs
- » Post-workout recovery
- » Mild athletic activity

Fiber Content

- » 67% *drirelease* (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- » Non-latex

Color Swatches



1. Effect of Compression Stockings on Physiological Responses and Running Performance in Division III Collegiate Cross-Country Runners During a Maximal Treadmill Test. Journal of Strength and Conditioning Research. Volume 28, Number 6, June 2014. Rider, Coughlin, Hew-Butler, Goslin.



401CX99
15–20mmHg
Closed Toe Calf Size X-Large in Black



401CM00
15–20mmHg
Closed Toe Calf Size M in White

401 Calf Sizing Chart

Women’s Shoe	<7.5	7.5–9.5	10–12	12.5–15	≥15.5
Men’s Shoe	<6.5	6.5–8.5	9–11	11.5–14	14.5–17
Ankle Circumference	7”–9” (18–23cm)	7.5”–10” (19–25cm)	8”–10.5” (20–27cm)	8.5”–11.5” (21–29cm)	8.5”–11.5” (21–29cm)
Size	S	M	L	X	2

Product Code Example: 401CL00

Series	Style	Size	Color
401	C	L	00

Styles and Colors

Closed Toe = ●	White 00	Black 99
Calf (C)	●	●

PERFORMANCE SLEEVES 412V

for women and men

Compression Level

20–30mmHg

Features and Benefits

- » Fabric resists deterioration in salt, chlorine, or fresh water
- » Helps reduce lactic acid buildup¹
- » Quick drying off
- » High-visibility reflective transfer

Best Suited for

- » Running, triathlons, endurance sports, and gym workouts
- » Team sports
- » Occupational wear

Fiber Content

- » 75% Nylon, 25% Spandex
- » Non-latex

Color Swatches



1. Effect of Compression Stockings on Physiological Responses and Running Performance in Division III Collegiate Cross-Country Runners During a Maximal Treadmill Test. Journal of Strength and Conditioning Research. Volume 28, Number 6, June 2014. Rider, Coughlin, Hew-Butler, Goslin.



412VS58
20–30mmHg
Calf Sleeve Size S in Purple



412VM50
20–30mmHg
Calf Sleeve Size M in Blue



412VL99
20–30mmHg
Calf Sleeve Size L in Black



412VM00
20–30mmHg
Calf Sleeve Size M in White

412V Calf Sleeve Sizing Chart

Ankle Circumference	7"–8.5" (18–21.5cm)	8.5"–10" (21.5–25.5cm)	10"–11.5" (25.5–29cm)	11.5"–13" (29–33cm)
Calf Circumference	11"–15.5" (28–39cm)	12.5"–17.5" (31.5–45cm)	14"–20" (35.5–51cm)	16"–23" (40.5–58cm)
Size	S	M	L	X

Product Code Example: 412VM58

Series	Style	Size	Color
412	V	S	58

Styles and Colors

Calf Sleeve = O	White 00	Blue 50	Purple 58*	Black 99
Calf Sleeve (V)	O	O	O	O

* Purple not available in X-Large

MERINO 421 OUTDOOR SOCK

for women and men

Compression Level

15–20mmHg

Features and Benefits

- » Soft and comfortable
- » Natural properties of wool make these socks odor resistant
- » Mesh flex zones for comfort and ventilation
- » Cushioned, padded soles protect against blisters
- » No compression in heel and instep for ease of donning
- » Thermoregulating for year-round comfort


Best Suited for

- » Tired, heavy, restless, aching legs
- » Outdoor sports and activities
- » Post-workout recovery

Fiber Content

- » 67% Fine Merino Wool,
26% Nylon, 7% Spandex
- » Non-latex

Color Swatches





421CL31
15–20mmHg
Closed Toe Calf Size L in Olive



421CM12
15–20mmHg
Closed Toe Calf Size M in Charcoal

421 Calf Sizing Chart

Women’s Shoe	<7.5	7.5–9.5	10–12	12.5–15
Men’s Shoe	<6.5	6.5–8.5	9–11	11.5–14
Ankle Circumference	7”–9” (18–23cm)	7.5”–10” (19–25cm)	8”–10.5” (20–27cm)	8.5”–11.5” (21–29cm)
Size	S	M	L	X

Product Code Example: 421CS12

Series	Style	Size	Color
421	C	M	12

Styles and Colors

Closed Toe = ●	Charcoal 12	Olive 31
Calf (C)	●	●

EVERSOFT 160 DIABETIC SOCK

for women and men

Compression Level

8-15mmHg

Features and Benefits

- » Thick, padded soles, heels, and toes for comfort and protection
- » Flat, low friction toe seam
- » Soft, stretchable knit-in band prevents pinching
- » Compression from ankle to knee, compression-free foot

Best Suited for

- » Individuals with diabetes
- » Individuals with limited dexterity or lower tolerance for compression
- » Individuals with sensitive skin

Fiber Content

- » 67% drilease (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- » Non-latex

Color Swatches





160CS00
8–15mmHg
Closed Toe Calf Size S in White



160CL00
8–15mmHg
Closed Toe Calf Size L in White

160 Calf Sizing Chart

Women’s Shoe	<7.5	7.5–9.5	10–12	12.5–15
Men’s Shoe	<6.5	6.5–8.5	9–11	11.5–14
Ankle Circumference	7”–9” (18–23cm)	7.5”–10” (19–25cm)	8”–10.5” (20–27cm)	8.5”–11.5” (21–29cm)
Size	S	M	L	X

Product Code Example: 160CS00

Series	Style	Size	Color
160	C	S	00

Styles and Colors

Closed Toe = ● White
00

Calf (C)	●
----------	---



ACCESSORIES

Donning Aids



SIGVARIS Doff N' Donner®

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use. Latex free.

Size	Product Code
One Size	583D100



Scan here to see
how to use Sigvaris
DOFF N' DONNER
bit.ly/sigdnd



SIGVARIS

Doff N' Donner Cone

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner. The bottom of the base firmly attaches to most clean, flat surfaces via suction.

Size	Product Code
One Size	583D200



MELANY®

MELANY from Sigvaris makes both donning and doffing compression garments simple. Its innovative one-piece design with adjustable handles was made especially for users with limited range of motion and finger dexterity.

Size	Calf Circumference	Product code
X-Small	25–30 cm	579SXS
Small	30–38 cm	579SSM
Medium	38–46 cm	579SMD
Large	46–56 cm	579SLG



Scan here to
see how to use
the MELANY
goto.sigvaris.com/bbiMdz

Sim-Slide (For open toe Styles)



Scan here to see how
to use the **Sim-Slide®**
bit.ly/SIGsyt

Size	Product Code
Small	587SLIDESM
Medium	587SLIDEMD
Large	587SLIDELG
X-Large	587SLIDEXL

Shoe size				
Women	≤4–6	6.5–8	8.5–9.5	≥10
Men	≤3–4.5	5–6.5	7–8	≥8.5
Size	S	M	L	XL

Magnide On/Off (For closed toe Styles)



Scan here to see how to
use the **Magnide On/Off®**
bit.ly/SIGmofyt

Size	Product Code
Medium	587MAGMD
Large	587MAGLG
X-Large	587MAGXL

Shoe size			
Women	≤8.5	9–12	≥ 12.5
Men	≤7.5	8–11	≥ 11.5
Size	M	L	XL

* XL suitable for extra-wide feet

Gloves



Rubber Gloves*
Ridged - The practical and protective way to apply and remove your compression stockings.

Package: 12 per Pack	
Size	Product Code
Small	591R400S
Medium	591R400M
Large	591R400L
X-Large	591R400X
ASST'D 4x S/M/L	591R400A

Rubber Glove and Latex-Free Glove Sizing Chart
Measure Circumference of the Hand at the Palm

6.5"-7" (16.5-18 cm)	7.5"-8" (19-20.5 cm)	8.5"-9" (21.5-23 cm)	9.5"-10" (24-25.5 cm)
Small	Medium	Large	X-Large

* **Caution:** Rubber Gloves contain natural rubber latex that may cause allergic reaction



Latex-Free Gloves
Gripping donning gloves recommended for use in a latex-free environment.

Package: 12 per Pack	
Size	Product Code
Small	592R400S
Medium	592R400M
Large	592R400L
X-Large	592R400X

It Stays®



Roll-on body adhesive
A gentle roll-on body adhesive that washes off with water and may be used to hold compression garments in place against the skin. IT STAYS! is pliable and moves with the body.

Size	Product Code
One 2 oz Bottle (60ml)	5811300
Twelve 2 oz. Bottles (60ml)	5811350

Washing Solution



Washing Solution
Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for washing delicate garments!

Package Options	
Contains	Product Code
Six 16 oz. Bottles (473ml)	586W900
One 16 oz. Bottle (473ml)	586W950
Twelve 2 oz. Bottles (60ml)	586W975

Liners



Cotton Underliner
10-15 mmHg
Can be used as an underlayer to protect fragile skin, secure wound dressings, add 10-15 mmHg of compression, or create an easy glide surface to aid in putting on and taking off compression stockings.

Package: 2 per Pack	
Size	Product Code
Small / Short	588USS
Small / Long	588USL
Medium / Short	588UMS
Medium / Long	588UML
Large-XL / Short	588ULS
Large-XL / Long	588ULL

DEALER PORTAL

The SIGVARIS GROUP Dealer Portal is your online location for securely ordering SIGVARIS GROUP products and managing your account. Place and track your orders, review order history, contact Customer Care, and much more at sigvaris-online.com.

DEALER PORTAL BENEFITS

24/7 access-order any time!

Three ways to securely place an order online:
Catalog, Quick Entry or Excel Upload

Access to order history and tracking


Shipping options including drop ship and easy access to an address table for quick shipment of reorders

Check your account balance and print invoices

New product launch information

Current promotions

Inelastic compression order forms



Email sigvariscares@sigvaris.com and you will receive your active account within one business day!

HOW TO ORDER

Well Being products


Series Number	Product Style	Size	Color
120 = SHEER FASHION <i>(See individual pages for series number)</i>	N = Thigh-high <i>(See individual pages for available styles)</i>	A = Small <i>(See individual pages for sizing)</i>	33 = Honey <i>(See individual pages for available colors)</i>
120	N	A	33

To order **ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, MERINO OUTDOOR SOCKS, and EVERSOFIT DIABETIC SOCKS:**

Size is only one letter and the letter for Gender is not used in the ordering code; however, to order SHEER FASHION in Open Toe styles, a letter O after the size is required in the ordering code.

Example: 421CL31 = MERINO OUTDOOR SOCK, 15–20 mmHg, Closed toe Calf, Size Large, Color Olive

SIGVARIS
GROUP



Stay connected with
SIGVARIS GROUP

Find a retailer:
Visit sigvaris.com

Become a dealer:
1.800.322.7744
sigvariscares@sigvaris.com

SIGVARIS, INC.
1119 Highway 74 South
Peachtree City, GA 30269 USA



SIGVARIS GROUP USA
sigvaris.com

Sigvaris, SIGVARIS GROUP, and MELANY are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, and in many countries worldwide. © Copyright 2023 SIGVARIS, Inc. Supima is a registered trademark of the Supima Association of America Corporation. All rights reserved. *drirelease* is Copyright © 2023 *drirelease*. *drirelease* is an Optimizer Brands company. All Rights Reserved.
LIT-CATALOG-WB 01.23